

ZX trolley

Handling cylinders has never been easier

Following our continuing innovations and improvements in cylinder design, we have now turned our attention towards improving the cylinder trolley. The ZX trolley from BOC Healthcare has a whole range of innovative features, which will ease the management of your medical gases and provide your patients with benefits as well.

Features and benefits

- Ergonomic, contemporary design in line with modern healthcare.
- Innovative retaining strap to secure the cylinder safely.
- Big anti-static wheels ensure easy and safe movement around the hospital.
- Quiet operation due to the unique flexible retaining strap – so that patients are not disturbed.
- Convenient location device for ENTONOX® demand valve or face mask.
- Bumpers and rubber base to help prevent damage and paint chipping.

The design details will benefit portering staff, clinical staff and patients. Research suggests that patients feel more comfortable within a ward with a pleasant, smart clinical environment. Therefore, well-designed wards and the equipment to match may have an influence on patient recovery time.



Assured quality

The trolley is designed by BOC (Medical Community Registered Design No. 000283437-001).

Spare parts available

Part no.	Part name
19310980	ZX cylinder trolley
19310995	Kit rubber strap assembly (includes handle, strap, screws and plates)
19310996	Kit wheels and fixing
19310997	Kit rubber base plate
19310998	Kit bumpers, bungs and screws
19311712	Safety strap for handle
19310982	Adhesive instruction label

This trolley comes with a comprehensive instruction manual inside the transport box, as well as an abbreviated user instruction label fixed to the trolley.

Specifications

Part	Material/specification
Frame	Steel
Handle	Glass filled nylon
Retaining strap	Neoprene
Dimensions	Height: 1075 mm Width: 330 mm Depth: 320 mm
Usage	Designed for indoor use
Cylinder weight	9.8 kg
Maximum load	20 kg (full cylinder weight) <small>Note: this recommendation is based on health and safety issues not trolley construction</small>
Finish	Powder coated
Bumpers	Polyethylene
Wheels	Anti-static 200mm solid rubber
Suitable cylinder	Only use cylinders with the following dimensions:
Sizes	Height: 930 – 940 mm Diameter: 140 – 143 mm Corresponds with the following BOC cylinder sizes: F, HX, EX and ZX
Cleaning	Soap and water with soft cloth or brush
Disinfecting	Virkon solution

To find out more or place your order today contact:

08456 447 764

or for 24 hour support visit:

www.bochealthcare.co.uk



BOC Healthcare, Great Britain & N. Ireland
Customer Service Centre, Priestley Road,
Worsley, Manchester M28 2UT, United Kingdom
Tel 08456 447 764, Fax 0845 087 0062
bochealthcare-uk@boc.com, www.bochealthcare.co.uk

BOC Healthcare, Republic of Ireland
PO Box 201, Bluebell, Dublin 12, Republic of Ireland
Tel 01 409 1827, Fax 01 409 1854
healthcareinfo.ie@boc.com, www.BOCOnline.ie