

Emergency pain relief

Small, lightweight cylinder

ENTONOX® EA cylinders give you more gas than any other cylinder of the same size.

The new EA cylinder has the same features and benefits of all our tried and trusted integral valve cylinders:

- integrated pressure regulator with no need to waste time fitting separate regulators
- easy on/off hand wheel operation for improved safety and time saving
- “live” contents gauge, indicating contents level even when the cylinder is turned off
- 350 litres are contained in the ENTONOX® EA cylinder giving approximately 58 minutes of pain relief¹
- ideal in first response situations to provide pain relief until the patient can be transported to the ambulance.

¹ based on an average adult consumption rate of 6 litres per minute.

Cylinder specifications

Cylinder order code	211-EA
Nominal contents (litres)	350
Nominal cylinder pressure (bar)	217
Nominal outlet pressure (bar)	4
Valve outlet pressure connection	ENTONOX® BS 5682 Schrader
Valve operation	Handwheel
Flow-rate (litres/min)	Schrader: 40
Dimensions (inc. valve)	366 x 85
L x D (mm)	
Water capacity (litres)	1.0
Nominal weight full (kg)	2.4

To find out more contact
0800 111 333 24 hours a day, or visit
www.bohealthcare.co.uk



BOC Healthcare, Great Britain & N. Ireland
Customer Service Centre, Priestley Road,
Worsley, Manchester M28 2UT, United Kingdom
Tel 0800 111 333, Fax 0800 111 555
bohealthcare-uk@boc.com, www.bohealthcare.co.uk

BOC Healthcare, Republic of Ireland
PO Box 201, Bluebell, Dublin 12, Republic of Ireland
Tel 01 409 1827, Fax 01 409 1854
healthcareinfo.ie@boc.com, www.BOCOnline.ie